


IOB WARMING UNIT

The Best Choice for Body Temperature Management to Combat Hypothermia

FDA CE

**SAFE
RELIABLE
EASY TO USE**

Rolling Stand Available



IOB MEDICAL

Copyright ©2014-2018, IOB Medical Inc. All Right Reserved.



Specifications

IOB Warming Unit IOB 505



Dimensions: H28cm x D22cm x W22cm

Weight: 5.4 kg

Relative Noise Level: <55 dBA

Operating Temperature:

High: $43^{\circ} \pm 2^{\circ}\text{C}$ ($109.4^{\circ} \pm 5.4^{\circ}\text{F}$)

Med: $38^{\circ} \pm 2^{\circ}\text{C}$ ($100.4^{\circ} \pm 5.4^{\circ}\text{F}$)

Low: $32^{\circ} \pm 2^{\circ}\text{C}$ ($89.6^{\circ} \pm 5.4^{\circ}\text{F}$)

Device Rating:

110-120 VAC, 60 Hz, 10 A

220-240 VAC, 50 Hz, 8 A

Filter: 0.2 μm High-efficiency Filter, 500 hours

Standards: IEC 60601-1, IEC 60601-1-2, IEC 60601-1-8, IEC 80601-2-35




Intraoperative Warming Blankets

Upper Body

Model: 003


195cm x 64cm



Lower Body

Model: IOB-002

142cm x 120 cm



Torso

Model: IOB-001

142cm x 120cm



Full Body Surgical

Model: IOB-008

210cm x 120cm



Underbody Warming Blankets

Pediatric Underbody

Model: IOB-005

91cm x 80cm



Large Pediatric Underbody

Model: IOB-009

160cm x 80cm



Spinal Underbody

Model: IOB-010


215cm x 100cm



Lithotomy Underbody

Model: IOB-011

170cm x 100cm



Full Access Underbody

Model: IOB-006

215cm x 100cm




Pediatric Warming Blankets

Pediatric Lower Body

Model: IOB-012

142cm x 100cm



Pediatric Long

Model: IOB-014

110 x 64 cm

(200cm x 17cm)



Pediatric Full Body

Model: IOB-007

170cm x 100cm



Pediatric Underbody

Model: IOB-005

91cm x 80cm



Large Pediatric Underbody

Model: IOB-009

160cm x 80cm



Postoperative Warming Blankets

Adult Full Body

Model: IOB-004

210cm x 120cm



Pediatric Full Body

Model: IOB-007

170cm x 100cm



Other Warming Blankets

Cath Lab

Model: IOB-015

180 x 74cm

(380cm x 17cm)



Surgical Access

Model: IOB-016

210cm x 120cm



Cardiac Access

Model - 020

142cm-x 120cm

