

PRAGATHI MEDICAL SYSTEMS

OXYGEN HOOD



FEATURES

- Single Piece, Round Shape, Unbreakable, Autoclavable polycarbonate material.
- Trauma free silicon neck adjustment strap
- Bilateral oxygen nozzle (prevents direct flow of cold oxygen to patients head)
- Width: 8.5" ; Height: 7.5" (medium), 6.5" (small)

E-mail: sales@pragathineo-natal.com

Phone: +91-9844111792
+91-9449764764



13485:2016