

PRAGATHI MEDICAL SYSTEMS



Baby Weighing Scale 3 in 1 BWS 2121



FEATURES

- Sturdy Platform to weigh an adult
- High bright red 2 displays for strain free reading
- Zeroing facility
- Operating Voltage: 230v-50Hz
- Dimensions:
 - Base - 45 cm x 55 cm
 - Height - 112 cm
- Battery backup - 5 hours.

Capacity	Resolution	Accuracy
200 Kg	50 Gm	50 Gm

E-mail: sales@pragathineo-natal.com
Phone: +91-9844111792
+91-9449764764



PRAGATHI MEDICAL SYSTEMS



Baby Weighing Scale 2 in 1 BWS 2122

FEATURES

- Sturdy Platform to weigh an adult
- High bright red 2 displays for strain free reading
- Zeroing facility
- Operating Voltage: 230v-50Hz
- Dimensions:
 - Base - 45 cm x 55 cm
 - Height - 112 cm
- Battery backup - 5 hours.



Capacity	Resolution	Accuracy
200 Kg	50 Gm	50 Gm

E-mail: sales@pragathineo-natal.com
Phone : +91-9844111792
+91-9449764764



BWS 2123

Baby Weighing Scale



Front View



Back View

Capacity	Resolution	Accuracy
30 Kgs	5 gm	± 5 gm

Features

- * Durable fibre baby tray
- * High bright display
Front & back
- * Operating voltage
230V-50HzAC
- * On & Off Key
- * Tare Functions
- * Auto Power Saver
- * Maintenance free
- * Battery Backup 5hour
Rechargeable type

PRAGATHI MEDICAL SYSTEMS

No. 49/1, Beside Ullas Theatre, Industrial Suburb, Yeshwanthpura,
Bengaluru - 560022

Ph: 9844111792, 9449764764